

## ***Procedure for making SU-8-2050 masters in the BMMSL***

### **Clean and spin**

- 1) Turn spinner ON and center wafer on 3" chuck.
- 2) Engage vacuum using the foot pedal. Make sure wafer is centered by setting the spin speed at a low RPM value and observing the wafer. Adjust as necessary.
- 3) Set speed to 3500 RPM and spray wafer with acetone twice followed by isopropyl alcohol (IPA) twice. Let the wafer spin for a few seconds.
- 4) Place wafer on a 125 °C hotplate. Let bake 5 min. Remove wafer and place on a clean wipe.
- 5) Turn hotplate OFF then ON and let cool to 65 °C. Set hotplate temperature to 65 °C.
- 6) Center wafer on 3" chuck.
- 7) Turn dispensing unit ON and place a grey dispensing tip on the SU-8 syringe. Attach the syringe to the dispensing unit.
- 8) Set dispensing unit pressure to 5 PSI and time to 60 s.
- 9) Use the foot pedal to start dispenser. Starting in the center of the wafer, dispense SU-8 in a swirl outward until timer stops.
- 10) Disengage vacuum on the spinner. Set timer knob to 40 s.
- 11) Engage vacuum. Turn speed knob to 500 RPM and wait 10 s. SLOWLY (250 RPM s<sup>-1</sup>) increase spin speed up to 2000 RPM.
- 12) When wafer stops, turn spin speed down to 0 and timer up to 120. Disengage vacuum.
- 13) Engage vacuum, set the wafer spinning at a very slow speed (60 RPM), and remove edge bead in one of two ways: A) With a weigh spatula, or B) with an acetone-soaked swab.
- 14) Disengage vacuum, remove wafer, and wipe off any SU-8 tendrils or excess SU-8 on the bottom of the wafer.
- 15) Place wafer on a 4" wafer on the hotplate. Once all wafers are on the hotplate, set it to 95 °C for 35 min and cover with Al foil top. Enable AUTO OFF. It takes about 1 h for the hotplate to cool to room temperature.

### **Expose**

- 16) Turn nitrogen on in dressing room.
- 17) Once the hotplate is cool (step 16), carefully separate 3" wafers from 4" wafers using a razor blade.
- 18) On the exposure unit, build a sandwich with wafer, film, and filter glass. Make sure the emulsion (dull) side of the mask is down.
- 19) Set exposure energy to 1400 mJ/cm<sup>2</sup> (700 on the dial, two exposures each wafer).
- 20) Press EXPOSURE GAP then EXPOSE.
- 21) When the first exposure is done, press EXPOSURE GAP again, then MASK RAISE. Lift the mask frame slightly to reset the system.

- 22) Press EXPOSURE GAP then EXPOSE. When the second exposure is done, reset the system and remove your materials.
- 23) Don't forget to turn off the nitrogen tank in the dressing room!
- 24) Place wafer back on hotplate. When all wafers have been positioned, set hotplate to 95 °C for 45 min. Enable AUTO OFF.

## **Develop**

- 25)
- 26)

When done, separate 3'' wafers from 4'' wafers; pour just enough SU-8 developer in a crystallization dish to cover the wafer; place one wafer in dish and periodically swirl until all unexposed PR is removed; wait 3 additional minutes; on the second wafer, wait 5 min once all SU-8 is gone

Rinse with IPA; if white residue appears, rinse with acetone (2-3 sec) and then with IPA; blow dry with nitrogen

Bake on hotplate at 150C for 30 min

Last updated 1/20/2005  
Prepared for the Tiny Biotools Lab

## **Change log**